


How do I use my Metered- Dose Inhaler (MDI) with a spacer device?

CUE WORDS


1. **Remove** the caps.



2. **Shake** the MDI.



3. **Assemble** MDI mouthpiece in the spacer.



4. **Seal lips** around the SPACER DEVICE mouthpiece. **Breathe out.**


5. **Depress MDI** into spacer.

6. **Breathe in** slowly and deeply.

7. **Hold breathe** for 5 – 10 seconds.

8. **Breathe out.**

9. **Replace caps.**



THE  LUNG ASSOCIATION™
When you can't breathe, nothing else matters.

 and 
ASTHMA TRAINING & EDUCATOR COURSE COPD TRAINING & EDUCATOR COURSE

Metered dose inhaler with a spacer